


# FIV Ab/FELV Ag Rapid Diagnostic Test

The Biopanda FIV Ab + FELV Ag Combined Rapid Test qualitatively detects Feline Immunodeficiency Virus Antibody (FIV Ab) and Feline Leukaemia Virus Antigen (FELV Ag) from a drop of cat's whole blood, serum, or plasma specimen.



## Test Features

- Easy to use
- Rapid results in only 5-10 minutes
- Results are read visually.
- Storage: 2-30°C



## Ordering Information

PRODUCT NAME	CATALOGUE No	KIT SIZE
FIV Ab/FELV Ag Rapid Test	RAPG-FIF-001	10 Tests
Feline Immunodeficiency Virus Ab Rapid Test	RAPG-FIV-001	10 Tests
Feline Leukaemia Virus Ag Rapid Test	RAPG-FELV-001	10 Tests

Email: [info@biopanda.co.uk](mailto:info@biopanda.co.uk)  
 Web: [www.biopanda.co.uk](http://www.biopanda.co.uk)

### Biopanda Reagents Ltd.

Unit 14, Carrowreagh Business  
 Park Carrowreagh Road  
 Belfast, BT16 1QQ, United  
 Kingdom  
 Tel : +44 (0)28 95438774