



BIODX VETERINARY DIAGNOSTIC SERVICES

Postal Address: PO Box 105, Macarthur Square, Campbelltown NSW 2560
Street Address: Suite 3, 251 Queen Street, Campbelltown, NSW 2650
P: 0402152194; E: info@biodx.com.au; W: www.biodx.com.au

Lab Use Only

Lab ID:
Date received:

WORM TEST REQUEST FORM

Submitter

Sender/ Contact person:

Full Address:

*Phone: *Email (print Clearly):

Name and Address of Client/Farm Being Tested (if different to above):

How do you wish to receive results? Email WhatsApp (Number?)

Animal(s) Species: Horse Donkey Cattle Sheep Goat Alpaca Other:

Sample #	1	2	3	4	5	6
Animal/Group ID						
Age/Class						
Date of last Worming						
Name of Wormer product						
Is the Animal Pregnant						
Other Note						

Sample collection date:..... Collection time:.....

Test Requests (tick box)

WORM EGG COUNT: Equine Standard FEC Equine Standard FEC (BULK) Ruminant Standard FEC
 Ruminant Standard FEC + Coccidia Test Ruminant Standard FEC (BULK) FERT Coccidia Only

Sample Collection and Mailing Instructions - See overpage

Payment: Please fill out the form overpage

Credit Card Purchased through website - Provide order number

Test Kit Orders Please send me more sample test submission kits. No. required:

Shipping Address (if different to above):

Payment Form

Prices are available on the website (biodx.com.au/wormtest)

Amount paid or Authorised: \$..... Date: ---- / ---- / ----

Payment Method:

Credit Card

 VISA  MASTERCARD

Name on Card:

Number: _____ - _____ - _____ - _____

Expiry date: ____ / ____ CCV: ____ ____

Send me an tax invoice

Sample Collection

1. Faeces: Place at least one tablespoonful (5-10 grams) of fresh faeces per animal in individual small bag or container. If using container with caps do not fill to the top, or they will leak.
2. For bulk tests: In separate bags, collect golf ball-sized samples from maximum ~10 individuals for each group. We'll combine them before testing.

Sample Shipment

3. Mailing: Seal all bags carefully, preferably pack into a second sealed bag, and place in mailing satchel. Include the completed submission form(s) but take care to ensure they cannot be damaged by leaking faeces or moisture. Postage costs will be minimised if all samples are sent in one parcel, consider including an ice brick if the weather is hot. Try Express Post samples are priority processed. If you elect to use this service, simply place your samples inside an Express Post satchel. Retain the tracking number sticker so you can trace your samples in transit. If possible, try to submit directly to our lab which save time and increase accuracy.
4. **SHIPPING ADDRESS: BIODX, SUITE 3 (Level 1), 251 QUEEN STREET, CAMPBELLTOWN, NSW 2560**

Frequently Asked Questions

Should I request a bulk test, or individual counts?

Equine: Most equine samples are analysed as individual counts, which indicate the parasite load in each animal. Bulk tests can be considered for groups of animals kept together, as they are cheaper and give you a good idea of the overall parasite load in the group.

Ruminant: A bulk test is cheaper and will give you a good idea of the overall parasite load for your herd/mob. Individual counts give greater detail about the distribution of parasites in a mob (whether some animals are more infected than others) and are also recommended in situations where infected animal kept separately.

When should I re-test?

Equine: If your counts are high, you should re-test your animals 2+ weeks after treatment to confirm the drench has worked. If your worm egg counts are low, you should consider re-testing in approximately 8 weeks, or sooner after seasonal changes (especially increased rainfall) or times of increased animal stress. If you're concerned that worms may be increasing, then re-test sooner, after approximately 4 weeks.

Ruminant: If your counts are high, you should re-test your animals approximately 10-14 days after treatment. If your worm egg counts are low, you should consider re-testing in approximately 6-8 weeks, or sooner after seasonal changes (especially increased rainfall) or times of increased animal stress. If you're concerned that worms may be increasing, then re-test sooner, after approximately 4 weeks

What included in Prepaid test Kits and Price?

Prepaid kits include: 1. One Reply paid sachet 2. 10 Ziplock sample bags 3. One pair of gloves and 4. A submission form. Price - \$18 per kits (free postage).